# Instruction For Use

## PANTHER-S(Sterile Spinal Cage)

### 1. Instruction For Use

#### 1.1 Language Selection

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### ◆ INTERBODY FUSION DEVICE (STERILE SPINAL CAGE)

This product is interbody fusion devices for intervertebral body fusion procedures in skeletally mature patients with degenerative disc disease(DDD) of the lumbar spine at one or two contiguous levels from L2 to S1. The products are divided into a general Full Mesh Type and a Graft Hole Type that can be implanted by filling artificial bones.

### ◆ INDICATION FOR USE

This product used to treat structural abnormalities caused by degenerative intervertebral disc of lumbar spine. This product is manufactured using a 3D printer.

### ◆ PRECAUTIONS BEFORE USE

- 1) Do not damage the product.
- 2) It should be treated by a professional medical staff and should be followed by the correct method of use.
- 3) Before the procedure, the user checks the patient's condition and selects an appropriate implant.
- 4) This product is supplied with Gamma radiation (ISO11137-1,2,3), re-sterilization is not allowed, and The damage and expiration date of the product's packaging must be checked before surgery.
- 5) Check whether there is any problem with the product or its packaging condition and whether it is damaged or damaged.
- 6) Use other than for purposes of use is prohibited.

### ◆ METHOD OF USE

- 1) Posterior lumbar (PLIF), Anterior lumbar (ALIF), Transforaminal lumbar (TLIF), Olique lumbar(OLIF), The treatment method determines the shape of the implant. Use the X-ray and C-arm to determine the level of disc.
- 2) After incision, open the affected area (upper and lower body) using the opening device.
- 3) Use Shaver to remove damaged disc altogether or partially to make a space for the implant to be inserted.
- 4) Use Trial to check the height and angle of the space from which the disc has been removed and determine the appropriate size of the product.
- 5) Insert the implant using the Cage holder.
- 6) After final verification of the fixation and calibration status of the inserted implant, suture the surgical site.

### ◆ MANAGEMENT METHOD AFTER USE

This product is sterile disposable medical devices and is prohibited from reusing them.

### ◆ STORAGE CONDITION

Room temperature (1~30°C)

### ◆ WARNING

- 1) Only surgeons skilled in surgery can use it.
- 2) The product has a potential risk of requiring additional surgical treatment below.
  - ① Reoperation due to failure of bone fusion
  - ② dislocation
  - ③ Damage to the product
  - ④ a secondary fracture
  - ⑤ Infection
  - ⑥ damage to blood vessels / nerves / organs
- 3) This product is prohibited if the patient exhibits an allergic or foreign reaction to the raw material.
- 4) Check the expiration date of this product before using it.
- 5) Do not use products with damage to packaging or products.
- 6) This product is disposable medical devices and is prohibited from reusing them.
- 7) This product is supplied sterilized so re-sterilization is prohibited.
- 8) For overweight or obese patients, surgery should be considered depending on the patient's condition, as it can cause excessive stress on the product, causing failure of device fixation or damage to the device itself.

### ◆ CONTRAINDICATIONS

- Patients with active infection or spinal inflammation
- Patients with inflammatory degenerative diseases
- Hepatitis or HIV-infected patients
- In the event of severe burns, contusions, or inflammation to surrounding tissue
- Patients with allergic symptoms to metallic substances
- Infection fracture patients
- Overweight patients
- In case a pregnant woman has not confirmed the pregnancy
- Osteoporosis patient
- Those who cannot follow post-operative guidelines during post-operative recovery
- In case of abnormal loading or auxiliary reaction due to neurological or osteosarcomic diseases
- Bone immature patients
- Massive inflammation
- Multiple-sided infection due to hemorrhagic nature in the implant site.
- When the process of bone absorption or osteoporosis spreads at a rapid rate in radiological considerations.
- A person whose diagnostic name is outside of the service category and whose Concentrator determines that it cannot be use.

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The above contraindications can be relative or absolute and must be consulted when the procedure is determined by the surgeon. However, contraindications are not limited to this list.

### ◆ GENERAL CONDITIONS OF USE

- 1) This product must be treated only by a surgeon with specific training or experience.
- 2) Be careful of infection during surgery.
- 3) This product does not ensure exposure to the MR environment as safety and compatibility assessments, heating and movement are not tested in the magnetic resonance (MR) environment.
- 4) The user should discuss with the patient all physical and mental limitations of the use of this implant.
- 5) The user must provide information to the patient about that the implant may be broken or damaged by strenuous activity or trauma.
- 6) If the patient is involved in a job or activity that can exert undue pressure on the implant, the user must provide information to the patient about the potential risks and side effects of surgery.
- 7) The user must alert patients to the deterioration of synostosis due to smoking.

### ◆ SIDE EFFECT

- 1) Decrease in bone density due to changes in mechanical stress distribution shape
- 2) Damage to nerves or blood vessels if this product is moved or damaged due to improper activity
- 3) Damage to the product if the load is applied without sufficient bone fusion
- 4) Infection and inflammation of the epithelium or the heart
- 5) Neurological spinal cervical damage by trauma during surgery
- 6) Changes in spinal curvature and stiffness

### ◆ REPORTING SIDE EFFECT

In the event of side effects or accidents, it must be reported to the manufacturer.

### ◆ HANDLING

Correct handling of the implants is very important. The user should not scratch or crack the implant.

### ◆ COMPLAINTS

Upon receiving any complaints or complaints regarding the quality, identity, durability, reliability, validity, and/or performance of the product, the user should notify manufacturer or their local representative. Furthermore, if the instrument is malfunctioning or suspected to be malfunctioning, it should be immediately notified to manufacturer or its local representative. If products of manufacturer have been improperly operated or have caused or contributed to the patient's serious injury or death, contact its local representative or manufacturer by phone, fax, or document as soon as possible.

Please fill out the following when receiving all complaints.

- Manufacturing number of components (Lot, Number)
- Product name (Part name)
- Part, No
- Name and address to be contacted
- Explain the cause of the complaint in detail

For additional information or complaints, please contact:  
※Please check the attached label for product and manufacturer information.

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### ◆ Non-sterile spinal cage

This product is a porous form of cage manufactured with 3D printing technology using powder-based Select Laser Melting (SLM) method. A cage is transplanted between the spine of a patient with spinal disease to ensure sufficient fusion and mechanical stability. Depending on the shape and procedure of the product, it is classified as PLIF, TLIF, ALIF and OLIF, and it is classified as a typical Full Mesh Type and a Graft Hole Type, which can be implanted by filling in the artificial bone. Manufactured from Titanium Alloy (Ti-6Al-4V ELI Alloy), see below for details.

### ◆ Indication

Cage-type device used to treat structural abnormalities caused by degenerative intervertebral disc. It is a Interbody fusion cage manufactured using a 3D printer.

### ◆ PREPARATIONS BEFORE USE

- Specialized instruments must be used to ensure correct implantation.
- Instruments that have been overused or subjected to excessive force are more likely to be damaged when reused. Be sure to check the instrument for abrasion or damage before procedure.
- In addition, if corrosion, cracks, pinholes, etc. are found during use, they should be immediately discontinued and replaced with other instruments.
- In the event of a severe impact or drop on the floor, if damage or foreign substances are attached to the surgical instrument, it is impure and should be immediately discontinued and replaced with another instrument.

### ◆ STERILIZATION PROCEDURE RECOMMENDATION FOR NON-STERILE MEDICAL DEVICES INCLUDING IMPLANTS

Implants are provided in a non-sterile state, so sterilization by the user is required before use. In case of sterilization, sterilize it by placing it in a tray provided by manufacturer or in a sterile tray used by the hospital.

The following pressurized steam sterilization conditions are recommended by manufacturer;

- Only sterile products should exist in the operating room.
- To maintain aseptic warranty index of 10<sup>6</sup> (standard sterilization conditions), steam sterilization conditions used in hospitals are recommended in the

No	Sterilization method	Circular state	Temperature	Exposure time	Drying time
1	Steam sterilization	Pre-Vacuum	132°C (270°F)	4min	30min

following ways:

Sterilization in the paper filters, aseptic of the filter must be shall be verified before the pressure steam sterilization. If a user uses a different type of sterilization method, the user is responsible for it, and manufacturer assumes no responsibility.

### ◆ METHOD OF USE

The surgeon should be well informed of the surgical procedures, surgical instruments, and implant characteristics before surgery. Use the appropriate size and shape of the product for the patient, but do not overdo it for any other purpose.

- Posterior lumbar (PLIF), Anterior lumbar (ALIF), Transforaminal lumbar (TLIF), Olique lumbar(OLIF). The treatment method determines the shape of the implant and uses X-ray and C-arm to determine the level of disc to be operated on.
- After incision, open the affected area (upper and lower body) using the opening device.
- Use Shaver to remove damaged disks altogether or partially to free up space for implantation.
- Using the Trial, determine the correct size of the product by checking the height and angle between the vertebrae in the space where the disc has been removed.
- Insert the implant using the Cage holder.
- After final verification of the fixation and calibration status of the implanted implant,

### ◆ STORAGE AND MANAGEMENT METHOD AFTER USE

- Implants are disposable medical devices and are prohibited from reusing them.
- Implants are transported packaged and must not be damaged upon arrival.
- It must be stored in a clean, dry, and suitable temperature (Room temperature:1~30°C).

### ◆ PRE-OPERATIVE PRECAUTION

The implant of manufacturer should be used by a trained physician who is fully familiar with the required surgical methods. Doctors should take care not to apply improper stress to the spine or implants during the procedure. Follow the surgical procedures detailed in the surgical procedures provided by manufacturer. For example, when repositioning the instrument in place, avoid excessive force because it is likely to cause injury to the patient.

Care should be taken not to twist, hit or scratch the implant with the

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instrument to reduce the risk of damage. Extreme care is required when the instrument is used around organs, nerves, or vessels associated with life.

Special care is required when using the instrument in pediatrics.

Unless specifically stated on the label, instruments may be reused after decontamination, cleaning and sterilization.

### ◆ CONTRAINDICATIONS

Contraindications may be relative or absolute. The selection of certain implants must be carefully considered against the overall evaluation of the patient. Following particulars can reduce the possibility of their successful outcome.

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### ◆ GENERAL CONDITIONS OF USE

Transplantation of spinal cage is a procedure that requires specific skills due to the presence of a risk of serious injury to the patient and must be performed only by experienced surgeons, or trained, who are essential for the use of these systems. The information contained in the attached Product Manual is essential for the use of this implant, but not sufficient.

In other words, this information should not be considered a substitute for careful patient selection without considering the following:

- The skill or experience of a surgeon.
  - Professional judgment
  - Preoperative planning and selection of implants
  - Anatomical and biomechanical knowledge of the spine
  - Understanding the mechanical properties and raw materials of the implants used
  - Trainings and techniques of spinal surgery and other relevant implants for transplantation
  - Ensuring patients' cooperation in following appropriately specified postoperative care programs
  - Performing additional scheduled postoperative care
- In addition, this implant does not ensure exposure to the MR environment as safety and compatibility assessments, heating and movement are not tested in the magnetic resonance (MR) environment.

### ◆ INFORMATION FOR PATIENTS

The surgeon should discuss with the patient all physical and mental limitations of the use of this implant. This includes rehab, physical therapy, and wearing proper aids prescribed by the doctor. Specific discussions should discuss early lifting of heavy objects, the level of activity, and the need for regular treatment. The surgeon must provide information to the patient about the potential risks and side effects of surgery, and also warn them that the implant cannot reproduce the flexibility, strength, reliability, or durability of the actual normal healthy bone, and that the implant may be broken or damaged by strenuous activity or trauma, and may need to be replaced at a later date. If the patient is involved in a job or activity that can exert undue pressure on the implant (e.g., excessive walking, running, lifting heavy weights, or muscle strain), the surgeon must warn the patient before the procedure that the consequences of the procedure may result in failure. The surgeon must alert patients to the deterioration of bone mass due to smoking, and warn them of potential consequences. In patients with degenerative diseases, if there is a high level of disease at the time of transplantation, the life expectancy of implants can be significantly reduced.

### ◆ INSTRUMENTS

Specialized instruments are provided by manufacturer and must be used to ensure correct implantation. In rare cases, the instrument may be broken during the procedure. Instruments that have been overused or subjected to excessive force are more likely to be damaged when reused. Be sure to check the instrument for abrasion or damage before procedure.

### ◆ REUSE

Implants should never be reused. Although it may not seem to be damaged, the implants used may have defects and may have potential damage to their integrity that reduces their service life. It is recommended to verify (check) that

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the instrument is in good condition and working condition before being used for surgery.

### ◆ HANDLING

Correct handling of the implants is very important. The surgeon should not scratch or crack the implant.

### ◆ IMPLANT SELECTION AND USE

Choosing the right implant shape, size and design for each patient is critical to the success of the surgery. The doctor is responsible to do the right choice in accordance with the patient. Excessive weight of the patient may bring additional stress and strain to the implant, which can lead to metal fatigue and/or failure or deformation of the implant. The size and shape of the bone structure determine the shape, size and type of the implant. Once transplanted, the implant is affected by stress and deformation. This repetitive stress on implants is a consideration by the physician when selecting the implant, which applies not only during the procedure but also during the post-treatment period. In fact, stress and deformation on the implants can lead to metal fatigue loads and fractures of the implants before the bone grafts are fully adhesion. This may result in additional side effects or earlier removal of the implant.

Improper selection, placement, location selection, and fixation of implants can result in unusual pressure, which can reduce the life of the implant. Outlining or bending bars or plates is recommended only when absolutely necessary, depending on the surgical method of each system. Contouring of bars should only be done when appropriate tools are available. Incorrectly outlined bars and bars that have been over-coordinated several times should not be implanted. The surgeon should be well informed of the surgical procedures, surgical instruments, and implant characteristics before surgery. For additional procedural information, see manufacturer procedure Protocol. Regular and additional inspections are recommended to check the condition and location of the implant and the condition of the connecting bone.

### ◆ SYSTEM COMPATIBILITY

When a small amount of corrosion occurs in the implanted metal and alloys, contact with other metals can accelerate the corrosion process. The presence of corrosion can speed up the damage of the implant, and the amount of metal compounds that spread through the body system also increases. Rods, Hooks, Screws, Internal fixing instrument (implants) such as connectors must be made of the same type or suitable metal. Because each manufacturer uses different materials, different durability and manufacturing requirements, and different design variables, the components of the system should not be mixed with any other manufacturer's products. manufacturer not be liable for any adverse effects caused by the mixing of implants.

suture the procedure area.

### ◆ COMPLAINTS

Upon receiving any complaints or complaints regarding the quality, identity, durability, reliability, validity, and/or performance of the product, medical professionals should notify manufacturer or their local representative. Furthermore, if the instrument is malfunctioning or suspected to be malfunctioning, it should be immediately notified to manufacturer or its local representative. If products of manufacturer have been improperly operated or have caused or contributed to the patient's serious injury or death, contact its local representative or manufacturer by phone, fax, or document as soon as possible.

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-Manufacturing number of components (Lot. Number)

-Product name

-Name and address to be contacted

-Explain the cause of the complaint in detail

For additional information or complaints, please contact:

Manufacturer : Mantiz Co., LTD.  
Tel : 053-252-5833  
Fax : 053-252-5834

2022.03

MZSU-F003-005 (Rev.000)