# Instruction For Use

## 1. Instruction For Use

### 1.1 Language Selection

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Currently, this manual is available in Korean and English. Language selection is based on the overall content of countries with high usage and import potential. If a request or need to add a language is made, it will be reflected in the user manual.

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# Instruction For Use

## Interbody Fusion device (STERILE CERVICAL CAGE)

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### ◆ INTERBODY FUSION DEVICE (STERILE SPINAL CAGE)

This product is a sterilized radiation (Gamma) product. It is a porous cervical intervertebral fusion prosthesis manufactured by SLM (Selective Laser Melting) 3D printing technology using powder. It is used by transplanting between the intervertebral bodies to obtain optimal stability. It can be classified into Convex type, Lordotic type and Parallel type according to the shape of the product, and is subdivided into a general Full Mesh Type and a Graft Hole Type that can be implanted by filling an artificial bone. This product is manufactured from Titanium Alloy Powder (Ti-6Al-4V ELI Alloy), please refer to the following for details.

### ◆ INDICATION FOR USE

This product used to treat structural abnormalities caused by degenerative intervertebral disc of lumbar spine. This product is manufactured using a 3D printer.

### ◆ PRECAUTIONS BEFORE USE

- 1) Do not damage the product.
- 2) It should be treated by a professional medical staff and should be followed by the correct method of use.
- 3) Before the procedure, the user checks the patient's condition and selects an appropriate implant.
- 4) This product is supplied with Gamma radiation (ISO11137-1,2,3), re-sterilization is not allowed, and The damage and expiration date of the product's packaging must be checked before surgery.
- 5) Check whether there is any problem with the product or its packaging condition and whether it is damaged or damaged.
- 6) Use other than for purposes of use is prohibited.

### ◆ METHOD OF USE

- 1) Posterior lumbar (PLIF), Anterior lumbar (ALIF), Transforaminal lumbar (TLIF), Oblique lumbar(OLIF), The treatment method determines the shape of the implant. Use the X-ray and C-arm to determine the level of disc.
- 2) After incision, open the affected area (upper and lower body) using the opening device.
- 3) Use Shaver to remove damaged disc altogether or partially to make a space for the implant to be inserted.
- 4) Use Trial to check the height and angle of the space from which the disc has been removed and determine the appropriate size of the product.
- 5) Insert the implant using the Cage holder.
- 6) After final verification of the fixation and calibration status of the inserted implant, suture the surgical site.

### ◆ MANAGEMENT METHOD AFTER USE

This product is sterile disposable medical devices and is prohibited from reusing them.

### ◆ STORAGE CONDITION

Room temperature (1~30°C)

### ◆ WARNING

- 1) Only surgeons skilled in surgery can use it.
- 2) The product has a potential risk of requiring additional surgical treatment below.
  - ① Reoperation due to failure of bone fusion
  - ② dislocation
  - ③ Damage to the product
  - ④ a secondary fracture
  - ⑤ Infection
  - ⑥ damage to blood vessels / nerves / organs
- 3) This product is prohibited if the patient exhibits an allergic or foreign reaction to the raw material.
- 4) Check the expiration date of this product before using it.
- 5) Do not use products with damage to packaging or products.
- 6) This product is disposable medical devices and is prohibited from reusing them.
- 7) This product is supplied sterilized so re-sterilization is prohibited.
- 8) For overweight or obese patients, surgery should be considered depending on the patient's condition, as it can cause excessive stress on the product, causing failure of device fixation or damage to the device itself.

### ◆ CONTRAINDICATIONS

- Patients with active infection or spinal inflammation
- Patients with inflammatory degenerative diseases
- Hepatitis or HIV-infected patients
- In the event of severe burns, contusions, or inflammation to surrounding tissue
- Patients with allergic symptoms to metallic substances
- Infection fracture patients
- Overweight patients
- In case a pregnant woman has not confirmed the pregnancy
- Osteoporosis patient
- Those who cannot follow post-operative guidelines during post-operative recovery
- In case of abnormal loading or auxiliary reaction due to neurological or osteosarcomatous diseases

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- Bone immature patients
- Massive inflammation
- Multiple-sided infection due to hemorrhagic nature in the implant site.
- When the process of bone absorption or osteoporosis spreads at a rapid rate in radiological considerations.
- A person whose diagnostic name is outside of the service category and whose Concentrator determines that it cannot be use.

The above contraindications can be relative or absolute and must be consulted when the procedure is determined by the surgeon. However, contraindications are not limited to this list.

## ◆ GENERAL CONDITIONS OF USE

- 1) This product must be treated only by a surgeon with specific training or experience.
- 2) Be careful of infection during surgery.
- 3) This product does not ensure exposure to the MR environment as safety and compatibility assessments, heating and movement are not tested in the magnetic resonance (MR) environment.
- 4) The user should discuss with the patient all physical and mental limitations of the use of this implant.
- 5) The user must provide information to the patient about that the implant may be broken or damaged by strenuous activity or trauma.
- 6) If the patient is involved in a job or activity that can exert undue pressure on the implant, the user must provide information to the patient about the potential risks and side effects of surgery.
- 7) The user must alert patients to the deterioration of synostosis due to smoking.

## ◆ SIDE EFFECT

- 1) Decrease in bone density due to changes in mechanical stress distribution shape
- 2) Damage to nerves or blood vessels if this product is moved or damaged due to improper activity
- 3) Damage to the product if the load is applied without sufficient bone fusion
- 4) Infection and inflammation of the epithelium or the heart
- 5) Neurological spinal cervical damage by trauma during surgery
- 6) Changes in spinal curvature and stiffness

## ◆ REPORTING SIDE EFFECT

In the event of side effects or accidents, it must be reported to the manufacturer.

## ◆ HANDLING

Correct handling of the implants is very important. The user should not scratch or crack the implant.

## ◆ COMPLAINTS

Upon receiving any complaints or complaints regarding the quality, identity, durability, reliability, validity, and/or performance of the product, the user should notify manufacturer or their local representative. Furthermore, if the instrument is malfunctioning or suspected to be malfunctioning, it should be immediately notified to manufacturer or its local representative. If products of manufacturer have been improperly operated or have caused or contributed to the patient's serious injury or death, contact its local representative or manufacturer by phone, fax, or document as soon as possible.

Please fill out the following when receiving all complaints.

- Manufacturing number of components (Lot. Number)
- Product name (Part name)
- Part. No
- Name and address to be contacted
- Explain the cause of the complaint in detail

For additional information or complaints, please contact:  
※Please check the attached label for product and manufacturer information.

Manufacturer : Mantiz Co., LTD.  
Tel : +82 53-252-5833  
Fax : +82 53-252-5834