

# Instruction For Use

## STANPAD

### 1. Instruction For Use

#### 1.1 Language Selection

Currently, this manual is available in English. Language selection is based on the overall content of countries with high usage and import potential. If a request or need to add a language is made, it will be reflected in the user manual.

#### 1.1 Instruction For Use : English

# Instruction For Use

## STANPAD

### STERILE HEMOSTATIC GAUZE

#### ◆ STERILE HEMOSTATIC DRESSING

This product is a form of gauze containing kaolin and is used for the purpose of assisting hemostasis through manual compression in trauma and has a blue X-ray Detectable line. When a bleeding wound is compressed with this product, the kaolin in the dressing assimilates with small water molecules in the blood, leaving a combination of a high concentration of platelet factors and blood clotting factors in the wound area. So it helps to induce natural blood clotting mechanisms. It is used for bleeding management such as hemostasis by applying physical pressure through this principle. See below for details. See below for details.

#### ◆ PURPOSE OF USE

It is used for bleeding management such as hemostasis by applying physical pressure to local areas such as cuts, lacerations, and abrasions. It can also be used for the temporary treatment of surgical wounds and traumatic bleeding in dermatology after surgery and surgery.

#### ◆ PRECAUTIONS BEFORE USE

- This product is supplied with radiation sterilization (Gama radiation, ISO11137-1,2,3) so you must check the packaging for damage and validity of the product.
- Be sure to read the user manual carefully before using it.
- Observe all contraindications, warnings, and precautions listed in the user manual.

#### ◆ METHOD OF USE

- Unpack the product to take out hemostatic dressing.
- Cover or position hemostatic dressing completely in areas where bleeding occurs. All areas of bleeding must be in contact.
- Apply pressure until the bleeding stops. More than one dressing may be required to stop bleeding.
- If necessary, bandages can be applied over hemostatic dressing to maintain pressure.
- Dressing is removed smoothly by pulling the edges of the dressing within 24 hours. If it is difficult to remove the dressing, moisturize it with saline.

#### ◆ STORAGE AND MANAGEMENT METHOD AFTER USE

**This product is disposable and must not be reused or re-sterilized.**

#### ◆ WARNING

- If it suspect an allergic reaction, stop using it and take immediate medical action.
- This product is non-absorbable and should not be left in the wound.
- This product is a medical device and should not be used for other purposes than medical devices.

#### ◆ Contraindication

- Interaction Contraindication
  - When used in otolaryngology or obstetrics and gynecology, strong pressure is prohibited because it can cause swelling, ischemia, or necrosis of local tissue.
- Applicable Contraindication
  - This product prohibits the use of packaging that have already been opened or damaged.

#### ◆ Side Effect

- Infection

#### ◆ GENERAL PRECAUTIONS

- This product is disposable and must not be reused or re-sterilize.
- Reuse or re-sterilization may result in contamination of this product and/or risk of patient infection or cross-infection.
- Do not use the product if packaging is damaged.
- Do not use products that are past their expiration date

#### ◆ PRECAUTIONS OF APPLICATION

# Instruction For Use

## STANPAD

- Remove this product within 24 hours.
- If the product is cut to the desired size, the blue X-ray detectable line must be included in the area . Discard the rest parts.

### ◆ STORAGE AND HANDLING

- Storage method: room temperature storage (1~30°C)
- Expiration date: 8 years from the date of manufacture

### ◆ COMPLAINTS

Upon receiving any complaints or complaints regarding the quality, identity, durability, reliability, validity, and/or performance of the product, medical professionals should notify manufacturer or their local representative.

Please fill out the following when receiving all complaints.

- Manufacturing number of components (Lot. Number)
- Product name
- Name and address to be contacted
- Explain the cause of the complaint in detail

For additional information or complaints, please contact:

※Please check the attached label for product and manufacturer information.



**Manufacturer** : **Mantiz Co., LTD.**  
**Tel** : **053-252-5833**  
**Fax** : **053-252-5834**

2023.09

MZSU-F002-009 (Rev.005)